














ARTICULACIONES ACONSEJABLES

RECOMENDED ARTICULATIONS

اوامر يوصي بها

ARTICULATIONS CONSEILLÉES

ANTERIORES ANTERIORES سابق ANTÉRIEURES			POSTERIORES POSTERIORES وظيفة POSTÉRIEURES		
SUP. SUP. العالى SUP.	INFERIORES INFERIOR أدنى INFERIEURES		0°	20°	33°
62	4		29MZ	29	50
3M	5	4	29MZ	31	50
2D	2DD	4	29MZ	31	50
3M	5	4	29MZ	31	50
62	4		33M	32	52
3M	5	4	33M	33	52
62	4		33M	32	52
3M	5	4	33M	33	54
62	4		33M	33	54
3M	5	4	33M	33	54
62	4		33L	32	54
3M	5	4	33M	32	52
62	4		33M	32	52
3M	5	4	33L		34 / 54
62	4		33M	33	54
3M	5	4	33L	33	54
3M	5	4	33L		34 / 54
62	4				34 / 54
3M	5	4			34 / 54

En algunos casos se puede usar modelos especiales.

For short cases use specially moulds.

في الحالات الصغيرة ، استخدم قوالب خاصة.

Pour les petits cas, utilisez des moules spéciaux.



NORMAL		SHORT BITE	
0°	34 / 50 / 52 / 54	0°	29LL
20°	29 / 31 / 32 / 33	20°	F32
33°	29MZ / 33M / 33L		

LAS CIFRAS SIGNIFICAN: MEDIDAS (mm.)

THE FIGURES ARE MEANING: MEASURES (mm.)

(. الأرقام يعني : القياسات (مم)

LES NOMBRES SIGNIFIENT: MESURES (mm.)



62

W6: 40,2

LC: 8,5

WC: 7,6

Humedeciendo el talón de las piezas con monómero, se fortalece la unión con la base de la dentadura.

Moistening the heel of the pieces with monomer, the union with the base of the denture is strengthened.

مرهقة كعب القطع مع مونومر ، يتم تعزيز الاتحاد مع قاعدة طقم الأسنان.

En humidifiant le talon des pièces avec du monomère, l'union avec la base de la prothèse est renforcée.